

# COMPARISON OF AVAILABLE CONTINUOUS GLUCOSE MONITORS (CGM)



	DEXCOM G7 (10 DAY OR 15 DAY)	LIBRE 3+	INSTINCT	SIMPLERA	EVERSENSE
RECEIVER	Yes	Yes	No	No	No
PHONE APP	iPhone and Android	iPhone and Android	Medtronic app or app manager	Medtronic app	iPhone and Android – Must use
“REAL TIME” READINGS	Yes	Yes	Yes	No, must scan	Yes
HIGHS/LOWS ALARMS	Yes	Yes	Yes	Yes	Yes
INSULIN PUMP INTEGRATION	Ilet, Omnipod 5 (10 or 15 day), Tandem X2 or Mobi (10 day only)	Ilet, Tandem X2 or Mobi, Twist	Medtronic 780G	Medtronic 780G	Twist
WHERE SENSOR IS WORN	Arm	Arm	Arm	Arm	Implantable in arm
REPORTS AUTOMATICALLY SHARED WITH CLINIC	Yes if using phone app	Yes if using phone app	Yes if using Medtronic phone app	Yes if using Medtronic phone app	Yes
ABILITY TO SHARE RESULTS WITH FAMILY/FRIENDS	Yes if using phone app	Yes if using phone app	Yes if using Medtronic phone app	Yes if using Medtronic phone app	Yes
LENGTH OF USE	10 days or 15 days	15 days	15 days	6-7 days	Up to 1 year

\*Can use Dexcom Stelo and Abbott Lingo without a prescription. This information is for patient education only and does not represent endorsement of any vendors or products

**DIABETES CENTER OF EXCELLENCE**  
 AMBULATORY CARE CENTER (ACC), SECOND FLOOR  
 55 LAKE AVENUE NORTH, WORCESTER, MA 01655

NEW PATIENTS: **855-UMASS-MD** (855-862-7763)  
 EXISTING PATIENTS: **508-334-3206**

NS DIAB 0059 Rev. 04/14/26



**Diabetes Center of Excellence**  
[WWW.UMMHEALTH.ORG/DIABETES](http://WWW.UMMHEALTH.ORG/DIABETES)